

Orthopaedic Suction Cannula

A clear view of the product. A structured clinical reference.

AMT001

Product reference

Replaceable filter

Retains bone fragments while allowing fluid passage.



3 m tubing

Cannula, connectors and interchangeable tips.

DESIGNED FOR ORTHOPAEDIC SUCTION

A sterile, disposable suction set for drawing liquids during orthopaedic procedures where bone fragments are expected. The cannula can be disconnected from the tube, with replaceable filters and tips for the procedure.

Product specifications

Orthopaedic Suction Set / Reference AMT001

DIMENSIONS & SUPPLIED COMPONENTS	
Tube length	3 m
Tube inner diameter	7 mm
Cannula tip inner diameter	6 mm
Straight tip length	25 cm
Curved tip lengths	15 cm and 25 cm
Filter hole diameter	2.5 mm
Duration of use	Continuous use under 60 minutes

Materials

Cannula handle, filter and chamber: PC.
Tube, connector and cannula tip: PVC.

Packaging & product information

6 units per outer carton.
Ethylene oxide sterilisation.
ISO 13485 quality system. EU Class IIa.
Shelf life: 5 years.

Handling & storage

Store in a clean, dry place away from direct sunlight. Refer to the product instructions for use for clinical operation.

Product information transcribed from the supplied brochure. Portfolio layout refinement, October 2026.

